



## Run Info

Host Name	GXB02097 (localhost)
Position	X4
Experiment Name	VX06_H61211c
Sample ID	H61211c
Run ID	e8b913cf-61fa-4277-b24a-627b2f346601
Acquisition ID(s)	5960c71948de5ca9b4c9cb95a18998f9d0c61e4e, 4f5b739ce3229490970cb1331c2ab745f3388166
Flow Cell Id	FAL83583
Start Time	July 22, 10:54
Run Length	3d 0h 4m

## Run Summary

Reads Generated	12.68 M
Passed Bases	5.14 Gb
Failed Bases	2.86 Gb
Estimated Bases	8.01 Gb


## Run Parameters

Flow Cell Type	FLO-MIN106
Kit	SQK-DCS109
Initial bias voltage	-180 mV
FAST5 output	Enabled
FASTQ output	Enabled
BAM output	Disabled
Bulk file output	Disabled
Active channel selection	Enabled
Basecalling	Enabled
Specified run length	72 hours
FAST5 reads per file	4000
FAST5 output options	vbz_compress,fastq,raw
FASTQ reads per file	4000
FASTQ output options	compress
Mux scan period	1 hour 30 minutes
Reserved pores	0 %
Basecall model	High-accuracy basecalling
Read filtering	min_qscore=9


## Versions

MinKNOW	21.05.12
MinKNOW Core	4.3.7
Bream	6.2.5
Guppy	5.0.12

Cumulative Output Reads




Cumulative Output Bases




**Read Length Histogram Estimated Bases - Outliers Discarded**

Estimated N50: 745 b




**Read Length Histogram Basecalled Bases - Outliers Discarded**

Estimated N50: 770 b




**Read Length Histogram Estimated Bases**

Estimated N50: 745 b




**Read Length Histogram Basecalled Bases**


Estimated N50: 770 b




Duty Time Grouped




Duty time Categorised




Mux Scan Grouped




Mux Scan Categorised




Temperature History




Bias Voltage History



Translocation Speed




QScore





Disk Write Performance



## **Run Debug Messages**

- The sequencing run has finished, but basecalling may continue July 25, 10:58
- Mux scan for flow cell FAL83583 has found a total of 24 pores. 24 pores available for immediate sequencing July 25, 09:58
- Performing Mux Scan July 25, 09:56
- Mux scan for flow cell FAL83583 has found a total of 19 pores. 18 pores available for immediate sequencing July 25, 08:26
- Performing Mux Scan July 25, 08:23
- Mux scan for flow cell FAL83583 has found a total of 23 pores. 22 pores available for immediate sequencing July 25, 06:53
- Performing Mux Scan July 25, 06:51
- Mux scan for flow cell FAL83583 has found a total of 23 pores. 22 pores available for immediate sequencing July 25, 05:21
- Performing Mux Scan July 25, 05:18
- Mux scan for flow cell FAL83583 has found a total of 23 pores. 23 pores available for immediate sequencing July 25, 03:48
- Performing Mux Scan July 25, 03:46
- Mux scan for flow cell FAL83583 has found a total of 26 pores. 24 pores available for immediate sequencing July 25, 02:16
- Performing Mux Scan July 25, 02:13
- Mux scan for flow cell FAL83583 has found a total of 34 pores. 34 pores available for immediate sequencing July 25, 00:43
- Performing Mux Scan July 25, 00:41
- Mux scan for flow cell FAL83583 has found a total of 36 pores. 33 pores available for immediate sequencing July 24, 23:11
- Performing Mux Scan July 24, 23:08
- Mux scan for flow cell FAL83583 has found a total of 21 pores. 21 pores available for immediate sequencing July 24, 21:38
- Performing Mux Scan July 24, 21:36
- Mux scan for flow cell FAL83583 has found a total of 36 pores. 35 pores available for immediate sequencing July 24, 20:06
- Performing Mux Scan July 24, 20:03
- Mux scan for flow cell FAL83583 has found a total of 42 pores. 40 pores available for immediate sequencing July 24, 18:33
- Performing Mux Scan July 24, 18:31
- Mux scan for flow cell FAL83583 has found a total of 54 pores. 49 pores available for immediate sequencing July 24, 17:01
- Performing Mux Scan July 24, 16:58
- Mux scan for flow cell FAL83583 has found a total of 62 pores. 56 pores available for immediate sequencing July 24, 15:28
- Performing Mux Scan July 24, 15:26
- Mux scan for flow cell FAL83583 has found a total of 59 pores. 57 pores available for immediate sequencing July 24, 13:56
- Performing Mux Scan July 24, 13:53
- Mux scan for flow cell FAL83583 has found a total of 55 pores. 51 pores available for immediate sequencing July 24, 12:23
- Performing Mux Scan July 24, 12:21
- Mux scan for flow cell FAL83583 has found a total of 51 pores. 48 pores available for immediate sequencing July 24, 10:51
- Performing Mux Scan July 24, 10:48
- Mux scan for flow cell FAL83583 has found a total of 48 pores. 46 pores available for immediate sequencing July 24, 10:48

- sequencing July 24, 09:18
- Performing Mux Scan July 24, 09:16
- Mux scan for flow cell FAL83583 has found a total of 47 pores. 43 pores available for immediate sequencing July 24, 07:46
- Performing Mux Scan July 24, 07:43
- Mux scan for flow cell FAL83583 has found a total of 67 pores. 62 pores available for immediate sequencing July 24, 06:13
- Performing Mux Scan July 24, 06:11
- Mux scan for flow cell FAL83583 has found a total of 80 pores. 77 pores available for immediate sequencing July 24, 04:40
- Performing Mux Scan July 24, 04:38
- Mux scan for flow cell FAL83583 has found a total of 90 pores. 81 pores available for immediate sequencing July 24, 03:08
- Performing Mux Scan July 24, 03:06
- Mux scan for flow cell FAL83583 has found a total of 84 pores. 78 pores available for immediate sequencing July 24, 01:35
- Performing Mux Scan July 24, 01:33
- Mux scan for flow cell FAL83583 has found a total of 92 pores. 85 pores available for immediate sequencing July 24, 00:03
- Performing Mux Scan July 24, 00:00
- Mux scan for flow cell FAL83583 has found a total of 108 pores. 92 pores available for immediate sequencing July 23, 22:30
- Performing Mux Scan July 23, 22:28
- Mux scan for flow cell FAL83583 has found a total of 124 pores. 112 pores available for immediate sequencing July 23, 20:58
- Performing Mux Scan July 23, 20:55
- Mux scan for flow cell FAL83583 has found a total of 152 pores. 122 pores available for immediate sequencing July 23, 19:25
- Performing Mux Scan July 23, 19:23
- Mux scan for flow cell FAL83583 has found a total of 164 pores. 138 pores available for immediate sequencing July 23, 17:53
- Performing Mux Scan July 23, 17:50
- Mux scan for flow cell FAL83583 has found a total of 200 pores. 154 pores available for immediate sequencing July 23, 16:20
- Performing Mux Scan July 23, 16:18
- Mux scan for flow cell FAL83583 has found a total of 250 pores. 187 pores available for immediate sequencing July 23, 14:48
- Performing Mux Scan July 23, 14:45
- Mux scan for flow cell FAL83583 has found a total of 296 pores. 203 pores available for immediate sequencing July 23, 13:15
- Performing Mux Scan July 23, 13:13
- Mux scan for flow cell FAL83583 has found a total of 408 pores. 271 pores available for immediate sequencing July 23, 11:43
- Performing Mux Scan July 23, 11:40
- Mux scan for flow cell FAL83583 has found a total of 585 pores. 347 pores available for immediate sequencing July 23, 10:10
- Performing Mux Scan July 23, 10:08
- Mux scan for flow cell FAL83583 has found a total of 786 pores. 409 pores available for immediate sequencing July 23, 08:37
- Performing Mux Scan July 23, 08:35
- Mux scan for flow cell FAL83583 has found a total of 894 pores. 431 pores available for immediate sequencing July 23, 07:05

- Performing Mux Scan July 23, 07:02
- Mux scan for flow cell FAL83583 has found a total of 957 pores. 444 pores available for immediate sequencing July 23, 05:32
- Performing Mux Scan July 23, 05:30
- Mux scan for flow cell FAL83583 has found a total of 1031 pores. 457 pores available for immediate sequencing July 23, 03:59
- Performing Mux Scan July 23, 03:57
- Mux scan for flow cell FAL83583 has found a total of 1060 pores. 468 pores available for immediate sequencing July 23, 02:27
- Performing Mux Scan July 23, 02:24
- Mux scan for flow cell FAL83583 has found a total of 1099 pores. 471 pores available for immediate sequencing July 23, 00:54
- Performing Mux Scan July 23, 00:52
- Mux scan for flow cell FAL83583 has found a total of 1127 pores. 471 pores available for immediate sequencing July 22, 23:22
- Performing Mux Scan July 22, 23:19
- Mux scan for flow cell FAL83583 has found a total of 1180 pores. 481 pores available for immediate sequencing July 22, 21:49
- Performing Mux Scan July 22, 21:46
- Mux scan for flow cell FAL83583 has found a total of 1205 pores. 485 pores available for immediate sequencing July 22, 20:16
- Performing Mux Scan July 22, 20:14
- Mux scan for flow cell FAL83583 has found a total of 1242 pores. 487 pores available for immediate sequencing July 22, 18:44
- Performing Mux Scan July 22, 18:41
- Mux scan for flow cell FAL83583 has found a total of 1290 pores. 492 pores available for immediate sequencing July 22, 17:11
- Performing Mux Scan July 22, 17:09
- Mux scan for flow cell FAL83583 has found a total of 1333 pores. 491 pores available for immediate sequencing July 22, 15:38
- Performing Mux Scan July 22, 15:36
- Mux scan for flow cell FAL83583 has found a total of 1312 pores. 495 pores available for immediate sequencing July 22, 14:06
- Performing Mux Scan July 22, 14:03
- Mux scan for flow cell FAL83583 has found a total of 1356 pores. 495 pores available for immediate sequencing July 22, 12:33
- Performing Mux Scan July 22, 12:31
- Mux scan for flow cell FAL83583 has found a total of 1382 pores. 502 pores available for immediate sequencing July 22, 11:00
- Performing Mux Scan July 22, 10:58
- Starting sequencing procedure July 22, 10:58
- Waiting up to 300 seconds for temperature to stabilise at 34.0°C July 22, 10:54